## Disposal and Refurbish/Refoil of an ECD Cell

If it is necessary to dispose of an ECD cell, contact:

Nuclear Radiation Development Corp. 2937 Alt. Blvd. North Grand Island, NY 14072

Tel: (716) 773-7634 Fax: (716) 773-7744

for disposal instructions and current fees.

In addition, report the ECD cell disposal to:

PerkinElmer Instruments LLC Radiation Safety Officer 710 Bridgeport Ave. Shelton CT 06484

and

Nuclear Material Safety and Safeguard U.S. Nuclear Regulatory Commission Washington, DC 20251

and/or

your state and local agency, if applicable.